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Sequence - T1w-3D				
Study	CDIP - PCID / GEOMETRIC PHANTOM v1.0			
Vendor	GE	Philips	Philips	Siemens
Field Strength	3.0T	3.0T	3.0T	3.0T
Model	Discovery	Ingenia	Achieva	Trio
Version	23	R5	3.2.1	17
Sequence Name	3D FAST SPGR	3D TFE	3D TFE	3D MP-RAGE
Imaging Options	IRP - Asset	Fast (Sense)	Fast (Sense)	iPat
Pulse Timing				
TE (ms)	min full (2.932)	shortest (3.3)	shortest (3.3)	2.98
TR (ms)	min (6.66)	shortest (7.3)	shortest (7.3)	2300
Flip Angle (°)	11	9	9	9
TI (ms)	400	945	945	900
Scan Range				
FOV (in-plane) (mm)	256 x 256	256 x 248	256 x 248	256 x 256
Slice thickness (mm)	1	1	1	1
Gap between slices (mm)	0	0	0	0
No. Slices	224	224	224	224
Acquisition				
Orientation	Sagittal	Sagittal	Sagittal	Sagittal
Matrix size	256 x 256	256 x 248	256 x 248	256 x 256
Voxel size [L/R x A/P x I/S] (mm)	1 x 1 x 1	1 x 1 x 1	1 x 1 x 1	1 x 1 x 1
NEX	1	1	1	1
Acceleration factor (Parallel factor*)	2	2	2	2
Fold-Over direction	AP	AP	AP	AP
Reconstruction				
Matrix size	256	256	256	256
Voxel size [L/R x A/P x I/S] (mm)	1 x 1 x 1	1 x 1 x 1	1 x 1 x 1	1 x 1 x 1
Other				
Fat Suppression	None	None	None	None
Bandwidth	31.25kHz	228 Hz/px	228 Hz/px	240 Hz/px
Coil Type				
Head	x	x	x	x
Channel	8-12 (HNS)	15 (Head and Neck)	8	12
Timing				
Prescan Time+	00:30	00:30	00:30	00:30
Scan Time	08:00	08:00	08:00	08:00
Total Time (min)	08:30	08:30	08:30	08:30
Comments:	Tested 24-JUN-14	Tested 23-OCT-14	Tested 11-JUN-14	Tested 20-MAY-14

Sequence - DUAL PD/T2				
Study	CDIP - PCID / GEOMETRIC PHANTOM v1.0			
Vendor	GE	Philips	Philips	Siemens
Field Strength	3.0T	3.0T	3.0T	3.0T
Model	Discovery	Ingenia	Achieva	Trio
Version	23	R5	3.2.1	16
Sequence Name	FSE	TSE	TSE	TSE
Imaging Options	EDR, Asset	Sense	Sense	iPAT
Pulse Timing				
TE (ms)	min full (8.68) / 85	13/100	13/100	10/91
TR (ms)	3000	3000	3000	3000
Flip Angle (°)	125	90	90	165
TI (ms)	-	-	-	-
Scan Range				
FOV (in-plane) (mm)	240 x 240	240 x 240	240 x 240	240 x 240
Slice thickness (mm)	3	3	3	3
Gap between slices (mm)	0	0	0	0
No. Slices	48	48	48	48
Acquisition				
Orientation	Oblique axial	Oblique axial	Oblique axial	Oblique axial
Matrix size	256 x 256	256 x 254	256 x 254	256 x 256
Voxel size [L/R x A/P x I/S] (mm)	0.94 x 0.94 x 3	0.94 x 0.94 x 3	0.94 x 0.94 x 3	0.94 x 0.94 x 3
NEX	1	1	1	1
Acceleration factor (Parallel factor*)	2	1.7	2	2
Fold-Over direction	RL	RL	RL	RL
Reconstruction				
Matrix size	256	256	256	256
Voxel size [L/R x A/P x I/S] (mm)	0.94 x 0.94 x 3	0.94 x 0.94 x 3	0.94 x 0.94 x 3	0.94 x 0.94 x 3
Other				
Fat Suppression	Fat Sat	Fat Sat	Fat Sat	Fat Sat
Bandwidth	19.23 kHz	211 Hz/px	222 Hz/px	181 Hz/px
Coil Type				
Head	x	x	x	x
Channel	8-12 (HNS)	15 (Head and Neck)	8	12
Timing				
Prescan Time	00:30	00:30	00:30	00:30
Scan Time	02:43	06:12	05:24	05:17
Total Time (min)	03:13	06:42	05:54	11:34
Comments:	Phase FOV 0.75	Set parameter "image filter" at Weak		

For Siemens and GE only: Acquisition must be done twice. First one from superior part and a second to cover inferior part. Refer to technical manual for details.

Sequence - fMRI				
Study	CDIP - PCID / GEOMETRIC PHANTOM v1.0			
Vendor	GE	Philips	Philips	Siemens
Field Strength	3.0T	3.0T	3.0T	3.0T
Model	Discovery	Ingenia	Achieva	Trio
Version	22	R5	3.2.3	17
Sequence Name	fMRI EPI	FFE EPI	FFE EPI	fMRI EPI
Imaging Options	EDR, Fast Filtered, "eyes open"	CLEAR, SENSE, "eyes open"	CLEAR, SENSE, "eyes open"	"eyes open", iPat
Pulse Timing				
TE (ms)	30	30	30	30
TR (ms)	2400	2110	2110	2110
Flip Angle (°)	70	70	70	70
TI (ms)	-	-	-	-
Scan Range				
FOV (in-plane) (mm)	224 x 224	224 x 224	224 x 224	224 x 224
Slice thickness (mm)	3.5	3.5	3.5	3.5
Gap between slices (mm)	0	0	0	0
No. Slices	40	40	40	40
Acquisition				
Orientation	Oblique axial	Oblique axial	Oblique axial	Oblique axial
Matrix size	64 x 64	64 x 63	64 x 64	64 x 64
Voxel size [L/R x A/P x I/S] (mm)	3.5 x 3.5 x 3.5	3.5 x 3.5 x 3.5	3.5 x 3.5 x 3.5	3.5 x 3.5 x 3.5
NEX	1	1	1	1
Acceleration factor (Parallel factor*)	2	2.1	2	2
Fold-Over direction	AP	AP	AP	AP
Reconstruction				
Matrix size	64	64	64	64
Voxel size [L/R x A/P x I/S] (mm)	3.5 x 3.5 x 3.5	3.5 x 3.5 x 3.5	3.5 x 3.5 x 3.5	3.5 x 3.5 x 3.5
Other				
Fat Suppression	Fat Sat	Fat Sat. SPIR	Fat Sat. SPIR	Fat Sat.
Bandwidth (Hz/Px)	n/a	48.6	2441	2442
Number of acquisitions	100	100	100	100
Coil Type				
Head	x	x	x	x
Channel	8-12 (HNS)	15 (Head and Neck)	8	12
Timing				
Prescan Time+	00:30	00:30	00:30	00:30
Scan Time	04:00	04:00	04:00	04:00
Total Time (min)	04:30	04:30	04:30	04:30
Comments:		Please set parameter "dynamic stabilization" at : enhanced		

Sequence - T1w-3D				
Study	CDIP – PCID v3.7			
Vendor	GE	Philips	Philips	Siemens
Field Strength	3.0T	3.0T	3.0T	3.0T
Model	Discovery	Ingenia	Achieva	Trio
Version	23	R5	3.2.1	17
Sequence Name	3D FAST SPGR	3D TFE	3D TFE	3D MP-RAGE
Imaging Options	IrP - Asset	Fast (Sense)	Fast (Sense)	iPat
Pulse Timing				
TE (ms)	min full (2.932)	shortest (3.3)	shortest (3.3)	2.98
TR (ms)	min (6.66)	shortest (7.3)	shortest (7.3)	2300
Flip Angle (°)	11	9	9	9
TI (ms)	400	945	945	900
Scan Range				
FOV (in-plane) (mm)	256 x 256	256 x 248	256 x 248	256 x 256
Slice thickness (mm)	1	1	1	1
Gap between slices (mm)	0	0	0	0
No. Slices	180	180	180	192
Acquisition				
Orientation	Sagittal	Sagittal	Sagittal	Sagittal
Matrix size	256 x 256	256 x 248	256 x 248	256 x 256
Voxel size [L/R x A/P x I/S] (mm)	1 x 1 x 1	1 x 1 x 1	1 x 1 x 1	1 x 1 x 1
NEX	1	1	1	1
Acceleration factor (Parallel factor*)	2	2	2	2
Fold-Over direction	AP	AP	AP	AP
Reconstruction				
Matrix size	256	256	256	256
Voxel size [L/R x A/P x I/S] (mm)	1 x 1 x 1	1 x 1 x 1	1 x 1 x 1	1 x 1 x 1
Other				
Fat Suppression	None	None	None	None
Bandwidth	31.25kHz	228 Hz/px	228 Hz/px	240 Hz/px
Coil Type				
Head	x	x	x	x
Channel	8-12 (HNS)	15 (Head and Neck)	8	12
Timing				
Prescan Time+	00:30	00:30	00:30	00:30
Scan Time	04:52	06:20	06:17	05:21

Total Time (min)	05:22	06:50	06:47	05:51
-------------------------	--------------	--------------	--------------	--------------

I. Paramètres du protocole humain

Sequence - DUAL PD/T2				
Study	CDIP - PCID v3.7			
Vendor	GE	Philips	Philips	Siemens
Field Strength	3.0T	3.0T	3.0T	3.0T
Model	Discovery	Ingenia	Achieva	Trio
Version	23	R5	3.2.1	16
Sequence Name	FSE	TSE	TSE	TSE
Imaging Options	EDR, Asset	Sense	Sense	iPAT
Pulse Timing				
TE (ms)	min full (8.68) / 85	13/100	13/100	10/91
TR (ms)	3000	3000	3000	3000
Flip Angle (°)	125	90	90	165
TI (ms)	-	-	-	-
Scan Range				
FOV (in-plane) (mm)	240 x 240	240 x 240	240 x 240	240 x 240
Slice thickness (mm)	3	3	3	3
Gap between slices (mm)	0	0	0	0
No. Slices	48	48	48	48
Acquisition				
Orientation	Oblique axial	Oblique axial	Oblique axial	Oblique axial
Matrix size	256 x 256	256 x 254	256 x 254	256 x 256
Voxel size [L/R x A/P x I/S] (mm)	0.94 x 0.94 x 3	0.94 x 0.94 x 3	0.94 x 0.94 x 3	0.94 x 0.94 x 3
NEX	1	1	1	1
Acceleration factor (Parallel factor*)	2	1.7	2	2
Fold-Over direction	RL	RL	RL	RL
Reconstruction				
Matrix size	256	256	256	256
Voxel size [L/R x A/P x I/S] (mm)	0.94 x 0.94 x 3	0.94 x 0.94 x 3	0.94 x 0.94 x 3	0.94 x 0.94 x 3
Other				
Fat Suppression	Fat Sat	Fat Sat	Fat Sat	Fat Sat
Bandwidth	19.23 kHz	211 Hz/px	222 Hz/px	181 Hz/px
Coil Type				
Head	x	x	x	x
Channel	8-12 (HNS)	15 (Head and Neck)	8	12
Timing				
Prescan Time	00:30	00:30	00:30	00:30
Scan Time	02:43	06:12	05:24	05:17
Total Time (min)	03:13	06:42	05:54	05:47
Comments:	Phase FOV 0.75	Set parameter "image filter" at		

		Weak		
--	--	-------------	--	--

Sequence - 2D FLAIR				
Study	CDIP – PCID v3.7			
Vendor	GE	Philips	Philips	Siemens
Field Strength	3.0T	3.0T	3.0T	3.0T
Model	Discovery	Ingenia	Achieva	Trio
Version	22	R5	3.2.3	17
Sequence Name	2D FLAIR	2D FLAIR	2D FLAIR	2D TDF
Imaging Options	EDR, Asset, IR	Fast (Sense)	Fast (Sense)	iPat
Pulse Timing				
TE (ms)	140	125	125	123
TR (ms)	9000	9000	9000	9000
Flip Angle (°)	125	150	150	165
TI (ms)	2250	2500	2500	2500
Scan Range				
FOV (in-plane) (mm)	240 x 240	240 x 210	240 x 210	240 x 240
Slice thickness (mm)	3	3	3	3
Gap between slices (mm)	0	0	0	0
No. Slices	48	48	48	48
Acquisition				
Orientation	Oblique axial	Oblique axial	Oblique axial	Oblique axial
Matrix size	256 x 256	256x224	256x222	256 x 256
Voxel size [L/R x A/P x I/S] (mm)	0.94 x 0.94 x 3	0.94 x 0.95 x 3	0.94 x 0.95 x 3	0.94 x 0.94 x 3
NEX	1	1	1	1
Acceleration factor (Parallel factor*)	1	2	2	2
Fold-Over direction	RL	RL	RL	RL
Reconstruction				
Matrix size	256	256	256	256
Voxel size [L/R x A/P x I/S] (mm)	0.94 x 0.94 x 3	0.94 x 0.94 x 3	0.94 x 0.94 x 3	0.94 x 0.94 x 3
Other				
Fat Suppression	None	None	None	None
Bandwidth	27.78 kHz	164 Hz/px	242 Hz/px	222 Hz/px
Coil Type				
Head	x	x	x	x
Channel	8-12 (HNS)	15 (Head and Neck)	8	12
Timing				
Prescan Time+	00:30	00:30	00:30	00:30
Scan Time	04:32	04:48	04:12	04:05
Total Time (min)	05:02	05:18	04:42	04:35
Comments:	Flip angle should be set and not left to the scanner to decide	Set "image filter" parameter at Weak		

Sequence- T2*				
Study	CDIP – PCID v3.7			
Vendor	GE	Philips	Philips	Siemens
Field Strength	3.0T	3.0T	3.0T	3.0T
Model	Discovery	Ingenia	Achieva	Trio
Version	22	R5	3.2.3	
Sequence Name	GRE	FFE	FFE	GRE
Imaging Options	Asset	Sense	Sense	iPat
Pulse Timing				
TE (ms)	20	21	21	20
TR (ms)	650	650	650	650
Flip Angle (°)	20	20	20	20
TI (ms)	-	-	-	-
Scan Range				
FOV (in-plane) (mm)	240 x 240	240 x 240	240 x 240	240 x 240
Slice thickness (mm)	3	3	3	3
Gap between slices (mm)	0	0	0	0
No. Slices	48	48	48	48
Acquisition				
Orientation	Oblique axial	Oblique axial	Oblique axial	Oblique axial
Matrix size	256 x 256	256x256	256x254	256 x 256
Voxel size [L/R x A/P x I/S] (mm)	0.94 x 0.94 x 3	0.94 x 0.94 x 3	0.94 x 0.94 x 3	0.94 x 0.94 x 3
NEX	1	1	1	1
Acceleration factor (Parallel factor*)	1	2	2	2
Fold-Over direction	RL	RL	RL	RL
Reconstruction				
Matrix size	256	256	256	256
Voxel size [L/R x A/P x I/S] (mm)	0.94 x 0.94 x 3	0.94 x 0.94 x 3	0.94 x 0.94 x 3	0.94 x 0.94 x 3
Other				
Fat Suppression	None	None	None	None
Bandwidth (Hz/Px)	19.23 kHz	217 Hz/px	217 Hz/px	200 Hz/px
Coil Type				
Head	x	x	x	x
Channel	8-12 (HNS)	15 (Head and Neck)	8	12
Timing				
Prescan Time	00:30	00:30	00:30	00:30
Scan Time	02:15	04:13	04:18	03:04
Total Time (min)	02:45	04:43	04:48	03:34
Comments:	CV act_te = 20000	Set parameter "image filter" at Weak		
	Phase FOV 0.75			

Sequence - fMRI				
Study	CDIP – PCID v3.7			
Vendor	GE	Philips	Philips	Siemens
Field Strength	3.0T	3.0T	3.0T	3.0T
Model	Discovery	Ingenia	Achieva	Trio
Version	22	R5	3.2.3	17
Sequence Name	fMRI EPI	FFE EPI	FFE EPI	fMRI EPI
Imaging Options	EDR, Fast Filtered, "eyes open"	CLEAR, SENSE, "eyes open"	CLEAR, SENSE, "eyes open"	"eyes open", iPat
Pulse Timing				
TE (ms)	30	30	30	30
TR (ms)	2400	2110	2110	2110
Flip Angle (°)	70	70	70	70
TI (ms)	-	-	-	-
Scan Range				
FOV (in-plane) (mm)	224 x 224	224 x 224	224 x 224	224 x 224
Slice thickness (mm)	3.5	3.5	3.5	3.5
Gap between slices (mm)	0	0	0	0
No. Slices	40	40	40	40
Acquisition				
Orientation	Oblique axial	Oblique axial	Oblique axial	Oblique axial
Matrix size	64 x 64	64 x 63	64 x 64	64 x 64
Voxel size [L/R x A/P x I/S] (mm)	3.5 x 3.5 x 3.5	3.5 x 3.5 x 3.5	3.5 x 3.5 x 3.5	3.5 x 3.5 x 3.5
NEX	1	1	1	1
Acceleration factor (Parallel factor*)	2	2.1	2	2
Fold-Over direction	AP	AP	AP	AP
Reconstruction				
Matrix size	64	64	64	64
Voxel size [L/R x A/P x I/S] (mm)	3.5 x 3.5 x 3.5	3.5 x 3.5 x 3.5	3.5 x 3.5 x 3.5	3.5 x 3.5 x 3.5
Other				
Fat Suppression	Fat Sat	Fat Sat. SPIR	Fat Sat. SPIR	Fat Sat.
Bandwidth (Hz/Px)	n/a	48.6	2441	2442
Number of acquisitions	250	300	300	300
Coil Type				
Head	x	x	x	x
Channel	8-12 (HNS)	15 (Head and Neck)	8	12
Timing				
Prescan Time+	00:30	00:30	00:30	00:30
Scan Time	10:00	10:39	10:39	10:41
Total Time (min)	10:30	11:09	11:09	11:11
Comments:		Please set parameter "dynamic stabilization" at : enhanced		

Sequence – DTI					
CDIP – PCID v3.7					
Study					
Vendor	GE	Philips	Philips	Siemens	Siemens
Field Strength	3.0T	3.0T	3.0T	3.0T	3.0T
Model	Discovery	Ingenia	Achieva	Trio	Trio
Version	22	R5	3.2.1	17	17
Sequence Name	DWI	DWI	DWI	DWI	DWI
Imaging Options	ASS	Sense	Sense	iPat	iPat
Pulse Timing					
TE (ms)	min (87.3)	(shortest) 107	(shortest) 100	96	96
TR (ms)	9000	(shortest) 9976	(shortest) 9931	9400	9400
Flip Angle (°)	90	90	90	90	90
TI (ms)	-	-	-	-	-
Scan Range					
FOV (in-plane) (mm)	256 x 256	256 x 256	256 x 256	256 x 256	256 x 256
Slice thickness (mm)	2	2	2	2	2
Gap between slices (mm)	0	0	0	0	0
No. Slices	70	70	70	70	70
Acquisition					
Orientation	Oblique axial	Oblique axial - With eyes out of the way	Oblique axial - With eyes out of the way	Oblique axial - With eyes out of the way	Oblique axial
Matrix size	128 x 128	128 x 126	128 x 128	128 x 128	128 x 128
Voxel size [L/R x A/P x I/S] (mm)	2 x 2 x 2	2 x 2 x 2	2 x 2 x 2	2 x 2 x 2	2 x 2 x 2
NEX	1	1	1	1	1
Acceleration factor (Parallel factor**)	2	2	2	2	2
Fold-Over direction	AP	AP	AP	AP	AP
Reconstruction					
Matrix size	128	128	128	128	128
Voxel size [L/R x A/P x I/S] (mm)	2x2x2	2 x 2 x 2	2 x 2 x 2	2 x 2 x 2	2 x 2 x 2

Diffusion					
b-value 1	0	0	0	0	0
b-value 2	1000	1000	1000	1000	-
Number of directions	30	32	32	30	1
Other					
Fat Suppression	FatSat	FatSat / SPIR	FatSat / SPIR	FatSat	FatSat
Bandwidth (Hz/Px)	n/a			2056 Hz/px	2056 Hz/px
Coil Type					
Head	x	x	x	x	x
Channel	8-12 (HNS)	15 (Head and Neck)	8	12	12
Timing					
Prescan Time+	00:30	00:30	00:30	00:30	00:30
Scan Time	05:06	05:49	05:47	05:20	00:38
Total Time (min)	05:36	06:19	06:17	05:50	01:08
Comments:		Set parameter "directional resolution" at : OPT 32	Total 3 sequences (1 x DWI 30 directions + 2 x DWI 1 Direction b=0)		
			Total Time: 08:06	B02 +03 - DIFFUSION Mode: scan trace	

GE: T2 images (b0) = 3. If possible, place one at the beginning, one in the middle and one at the end / Make sure there is not interpolation. To avoid an automatic interpolation of data, please set control variable (CV). There are two ways of setting this CV, one with rhmethod = 1 or rhminsize = 128; both of these achieve the same end result.

Sequence - Task fMRI												
Study	CIMA-Q Add-on											
Vendor	Philips	Philips	Philips	Philips	Philips	Philips	Philips	Siemens	Siemens	Siemens	Siemens	
Field Strength	3.0T	3.0T	3.0T	3.0T	3.0T	3.0T	3.0T	3.0T	3.0T	3.0T	3.0T	
Model	Ingenia	Ingenia	Ingenia	Achieva	Achieva	Achieva	Trio	Trio	Trio	Trio	Trio	
Version	R5	R5	R5	3.2.3	3.2.3	3.2.3	17	17	17	17	17	
Sequence Name	FFE EPI	FFE EPI	B0 Map	FFE EPI	FFE EPI	B0 Map	fMRI EPI	fMRI EPI	fMRI EPI	fMRI EPI	GRE field mapping	
Imaging Options												
Pulse Timing												
TE (ms)	25	25	4.6	25	25	4.6	25	25	25	4.92 / 7.38		
TR (ms)	2500	2500	475	2500	2500	475	2500	2500	2500	476		
Flip Angle (°)	90	90	60	90	90	60	90	90	90	60		
TI (ms)	-	-	-	-	-	-	-	-	-	-		
Scan Range												
FOV (in-plane) (mm)	240 x 240	240 x 240	240 x 240	240 x 240	240 x 240	240 x 240	222 x 222	222 x 222	222 x 222	222 x 222		
Slice thickness (mm)	3	3	3	3	3	3	3	3	3	3		
Gap between slices (mm)	0.3	0.3	0.3	0.3	0.3	0.3	0	0	0	0		
No. Slices	41	41	41	41	41	41	41	41	41	45		
Acquisition												
Orientation	AC-PA minus 20°	AC-PA minus 20°	AC-PA minus 20°	AC-PA minus 20°	AC-PA minus 20°	AC-PA minus 20°	AC-PA minus 20°	AC-PA minus 20°	AC-PA minus 20°	AC-PA minus 20°	AC-PA minus 20°	
Matrix size	80 x 80	80 x 80	80 x 80	80x79	80x79	80 x 80	74	74	74	74	74	
Voxel size [L/R x A/P x I/S] (mm)	3 x 3 x 3	3 x 3 x 3	3 x 3 x 3	3 x 3 x 3	3 x 3 x 3	3 x 3 x 3	3 x 3 x 3	3 x 3 x 3	3 x 3 x 3	3 x 3 x 3	3 x 3 x 3	
NEX	1	1	1	1	1	1	1	1	1	1	1	
Acceleration factor (Parallel factor**)	1.8	1.8	1	1.3	1.3	1	1	1	1	1	1	
Fold-Over direction	AP	First AP / Second PA	AP	AP	First AP / Second PA	AP	AP	AP	First AP / Second PA	AP	AP	
Reconstruction												
Matrix size	80	80	80	64	64	80	74	74	74	74	74	
Voxel size [L/R x A/P x I/S]	3 x 3 x 3	3 x 3 x 3	3 x 3 x 3	3 x 3 x 3	3 x 3 x 3	3 x 3 x 3	3 x 3 x 3	3 x 3 x 3	3 x 3 x 3	3 x 3 x 3	3 x 3 x 3	

(mm)													
Other													
Fat Suppression	Fat Sat. SPIR	Fat Sat. SPIR	None	Fat Sat. SPIR	None	Fat Sat. SPIR	None	Fat Sat.	Fat Sat.	None			
Bandwidth (Hz/Px)	39.1	39.1	294.8	25.2	25.2	294.8	2502	2502	268				
Number of acquisitions	300	4	-	4	-	-	310	4	-				
Coil Type													
Head	X	X	X	X	X	X	X	X	X				
Channel	15 (Head and Neck)	15 (Head and Neck)	15 (Head and Neck)	8	8	8	12	12	12				
Timing													
Prescan Time+	00:30	00:30	00:30	00:30	00:30	00:30	00:30	00:30	00:30				
Scan Time	12:35	00:17	01:17	00:17	01:17	01:17	13:03	00:17	01:17				
Total Time (min)	13:05	00:47	01:47	00:47	01:47	01:47	13:33	00:47	01:47				
	Tested 22-AVR-15	Tested 22-AVR-15	Tested 22-AVR-15	Tested 10-AVR-15	-	-	Tested 10-MAR-15	Tested 10-MAR-15	Tested 10-MAR-15				
Comments:	Core task sequence need to be synchronized with e-Prime paradigm	Need to be done twice : Once with AP phase direction and after PA phase direction. Sequences for geometric distortion correction	Please acquire B0 Field Map after those three EPI sequences	Core task sequence need to be synchronized with e-Prime paradigm	Need to be done twice : Once with AP phase direction and after PA phase direction. Sequences for distortion correction	Please acquire B0 Field Map after those three EPI sequences	Core task sequence need to be synchronized with e-Prime paradigm	Need to be done twice : Once with AP phase direction and after PA phase direction. Sequences for distortion correction	Please acquire Gre Field Map after those three EPI sequences				

SCAN INFORMATION FORM

MRI Dataset Name: _____ **Type of Data:** Living Human Human Phantom Geometric Phantom
Site Name: _____ **Total Duration of Session:** _____ h _____ min
MRI Operator Name: _____ **Scan Date: (DD-MON-YYYY)** _____ - _____ - _____

Following section: To be completed by **MR technologist**
 Use **CDIP** approved sequences for all **CIMA-Q** scans. Please review the scan for motion and artifacts. Make sure **3D-T1w** and **FLAIR** are well acquired before scanning following sequences. Re-acquire if necessary. Refer to **CIMA-Q Manuel Scan procedure** for detailed instructions.
LORIS SUPPORTS ONLY CLASSIC DICOM, DO NOT EXPORT DATA IN ANY OTHER FORMAT (i.e ENHANCED DICOM)

Sequences	Acquired ?			Nb of attempts	Comments (e.g. subject woke up, repeated series #, motion)
	Yes	Partial	No		
3D-T1w					
FLAIR					
Dual PD/T2					
T2*					
DTI					
DWI B0 AP					
DWI B0 PA					
Connectivity fMRI					
Activation task					
Bold AP					
Bold PA					
Field Map					

Keep this copy for your own record.

Annex 6: Phantom Technical Information

PLEASE READ AND UNDERSTAND BEFORE USING THE PHANTOM

IMPORTANT INFORMATION AND INSTRUCTIONS

STORE THIS DOCUMENT WITH THE PHANTOM

Description:

The phantom container is made of poly-carbonate and was sealed with commercial water-resistant epoxy. The phantom contains 129 2x4 DUPLO® bricks, 15 2x2 bricks and 4 2x4 half-bricks and filled with 7.5 litres of aqueous solution containing in total approximately 0.15 grams of manganese chloride and 21 grams of sodium chloride. Bottom part of the phantom also contains a plastic vial filled with 5ml of another solution containing approximately less than 0.1 milligram of manganese chloride and 14 milligram of sodium chloride.

- Manganese chloride is a toxic substance and affects the human respiratory system when ingested in significant quantities. Please review the enclosed material safety data sheets.
- Sodium chloride may irritate skin and eyes on contact. Please review the enclosed material safety data sheets.

Handling:

The phantom is heavy (weighing around 8 kilograms) and may cause injury to personnel or damage to property if mishandled or dropped. The phantom container and phantom contents may shatter if dropped. Please exercise caution when carrying or handling the phantom.

- A person handling the phantom should be able to comfortably carry and handle the 8 kilogram phantom.
- Always carry the phantom from its bottom, holding it stably with two hands.
- Make sure the phantom is clean and dry before carrying or handling the phantom.
- Do not carry or manipulate the phantom by holding the container lid.
- Do not apply excessive pressure or force to the top, bottom, and sides of the phantom or the plastic cap on the lid.
- Do not pull, twist, or tug at the phantom container lid.
- Inspect phantom for leaks before and after each handling

Storage:

- Keep phantom in controlled temperature environment, ideally between 15 to 25 degrees Celsius (59°F – 77°F).
- Please store phantom in upright position.
- Do not stack items on top of the phantom.

MATERIAL SAFETY DATA SHEET

SECTION 1 - CHEMICAL PRODUCT AND COMPANY IDENTIFICATION

MSDS Identification:

- Key: 88295
- Name: Manganese chloride 1.0M solution

Catalog Numbers:

BP541-1, BP541-100

Synonyms:

- Manganese(II) chloride in solution; Manganese dichloride in aqueous solution; Manganous chloride solution.

Company Identification:

Fisher Scientific
1 Reagent Lane
Fairlawn, NJ 07410

For information, call:

- 201-796-7100

Emergency Number:

- 201-796-7100

For CHEMTREC assistance, call:

- 800-424-9300

For International CHEMTREC assistance, call:

- 703-527-3887

SECTION 2 - COMPOSITION, INFORMATION ON INGREDIENTS

CAS #	Chemical Name	%	EINECS #
7732-18-5	Water	87.5	231-791-2
7773-01-5	Manganese chloride	12.5	231-869-6

Text for R-phrases: see Section 16

Hazard Symbols: XN

Risk Phrases: 22

SECTION 3 - HAZARDS IDENTIFICATION

EMERGENCY OVERVIEW

Appearance: light pink clear liquid. Warning! May be harmful if swallowed. May cause central nervous system effects. Causes eye, skin, and respiratory tract irritation.

Target Organs: Central nervous system, lungs, reproductive system.

POTENTIAL HEALTH EFFECTS

Eye: Causes eye irritation. May cause chemical conjunctivitis.

Skin: Causes skin irritation.

Ingestion: May cause gastrointestinal irritation with nausea, vomiting and diarrhea. May be harmful if swallowed. May cause central nervous system effects and/or neurological effects. In high doses, manganese may increase anemia by interfering with iron absorption.

Inhalation: Causes respiratory tract irritation. The lowest exposure concentration of manganese at which early effects on the CNS and the lungs may occur is still unknown. However, once neurological signs are present, they tend to continue and worsen after exposure ends.

Chronic: Chronic inhalation or ingestion may result in manganism characterized by neurological symptoms such as headache, apathy, and weakness of the legs, followed by psychosis and neurological symptoms similar to those of Parkinson's disease. May impair fertility. Other chronic effects from inhaling high amounts of manganese include an increased incidence of cough and bronchitis and susceptibility to infectious lung disease.

SECTION 4 - FIRST AID MEASURES

Eyes: Immediately flush eyes with plenty of water for at least 15 minutes, occasionally lifting the upper and lower eyelids. Get medical aid.

Skin: Flush skin with plenty of water for at least 15 minutes while removing contaminated clothing and shoes. Get medical aid if irritation develops or persists. Wash clothing before reuse.

Ingestion: Do not induce vomiting. If victim is conscious and alert, give 2-4 cupfuls of milk or water. Never give anything by mouth to an unconscious person. Get medical aid immediately.

Inhalation: Remove from exposure and move to fresh air immediately. If not breathing, give artificial respiration. If breathing is difficult, give oxygen. Get medical aid.

Notes to Physician: Treat symptomatically and supportively. Persons with impaired respiratory function or anemia may be at increased risk to the hazards associated with this substance.

SECTION 5 - FIRE FIGHTING MEASURES

General Information: As in any fire, wear a self-contained breathing apparatus in pressure-demand, MSHA/NIOSH (approved or equivalent), and full protective gear. During a fire, irritating and highly toxic gases may be generated by thermal decomposition or combustion.

Extinguishing Media: Use extinguishing media most appropriate for the surrounding fire. Use water spray, dry chemical, carbon dioxide, or appropriate foam.

Autoignition Temperature: Not applicable.

Flash Point: Not applicable.

Explosion Limits, lower: Not available.

Explosion Limits, upper: Not available.

NFPA Rating: (estimated) Health: 2; Flammability: 0; Instability: 0

SECTION 6 - ACCIDENTAL RELEASE MEASURES

General Information: Use proper personal protective equipment as indicated in Section 8.

Spills/Leaks: Clean up spills immediately, observing precautions in the Protective Equipment section. Absorb spill using an absorbent, non-combustible material such as earth, sand, or vermiculite. Do not use combustible materials such as sawdust. Provide ventilation.

SECTION 7 - HANDLING AND STORAGE

Handling: Wash thoroughly after handling. Remove contaminated clothing and wash before reuse. Use only in a well-ventilated area. Avoid contact with eyes, skin, and clothing. Do not breathe dust, vapor, mist, or gas. Keep container tightly closed. Do not ingest or inhale.

Storage: Keep container closed when not in use. Store in a tightly closed container. Store in a cool, dry, well-ventilated area away from incompatible substances.

SECTION 8 - EXPOSURE CONTROLS, PERSONAL PROTECTION

Engineering Controls: Use process enclosure, local exhaust ventilation, or other engineering controls to control airborne levels below recommended exposure limits. Facilities storing or utilizing this material should be equipped with an eyewash facility and a safety shower.

EXPOSURE LIMITS

Chemical Name	ACGIH	NIOSH	OSHA - Final PELs
Water	none listed	none listed	none listed
Manganese chloride	0.2 mg/m ³ TWA (asMn) (listed under ** no name**).	1 mg/m ³ TWA (asMn) (listed under ** no name**); 500 mg/m ³ LDH (as Mn) (listed under ** no name**).	5 mg/m ³ Ceiling (as Mn) (listed under ** no name**).

OSHA Vacated PELs:

Water: No OSHA Vacated PELs are listed for this chemical.

Manganese chloride: No OSHA Vacated PELs are listed for this chemical.

PERSONAL PROTECTIVE EQUIPMENT

Eyes: Wear appropriate protective eyeglasses or chemical safety goggles as described by OSHA's eye and face protection regulations in 29 CFR 1910.133 or European Standard EN166.

Skin: Wear appropriate protective gloves to prevent skin exposure.

Clothing: Wear appropriate protective clothing to prevent skin exposure.

Respirators: A respiratory protection program that meets OSHA's 29 CFR 1910.134 and ANSI Z88.2 requirements or European Standard EN 149 must be followed whenever workplace conditions warrant a respirator's use.

SECTION 9 - PHYSICAL AND CHEMICAL PROPERTIES

Physical State: Clear liquid

Color: light pink

Odor: odorless

pH: Not available.

Vapor Pressure: Not available.

Vapor Density: Not available.

Evaporation Rate: Not available.

Viscosity: Not available.

Boiling Point: Not available.

Freezing/Melting Point: Not available.

Decomposition Temperature: Not available.

Solubility in water: Not available.

Specific Gravity/Density: Not available.

Molecular Formula: Solution

Molecular Weight: 0

SECTION 10 - STABILITY AND REACTIVITY

Chemical Stability: Stable under normal temperatures and pressures.

Conditions to Avoid: Excess heat.

Incompatibilities with Other Materials: Strong reducing agents, zinc, potassium, hydrogen peroxide, sodium.

Hazardous Decomposition Products: Hydrogen chloride, irritating and toxic fumes and gases, oxides of manganese.

Hazardous Polymerization: Will not occur.

SECTION 11 - TOXICOLOGICAL INFORMATION

RTECS#:

- CAS# 7732-18-5: ZC0110000
- CAS# 7773-01-5: OO9625000

LD50/LC50:

- CAS# 7732-18-5: Oral, rat: LD50 =>90 mL/kg.
- CAS# 7773-01-5: Oral, mouse: LD50 = 1031 mg/kg; Oral, mouse: LD50 = 450 mg/kg; Oral, rat: LD50 = 250 mg/kg. Not available.

Carcinogenicity:

Water -

- Not listed by ACGIH, IARC, or NTP.

Manganese chloride -

- Not listed by ACGIH, IARC, or NTP.

Epidemiology:

The U.S. EPA stated that epidemiological studies of inorganic manganese compounds in humans indicate effects on the respiratory system at levels below 1 mg/m³.

Teratogenicity:

No data available.

Reproductive Effects:

Men exposed to manganese dusts showed a decrease in fertility.

Neurotoxicity:

Manganese is neurotoxic.

Mutagenicity:

No data available.

Other Studies:

No data available.

SECTION 12 - ECOLOGICAL INFORMATION**SECTION 13 - DISPOSAL CONSIDERATIONS**

Chemical waste generators must determine whether a discarded chemical is classified as a hazardous waste. US EPA guidelines for the classification determination are listed in 40 CFR Parts 261.3. Additionally, waste generators must consult state and local hazardous waste regulations to ensure complete and accurate classification.

RCRA P-Series: None listed.

RCRA U-Series: None listed.

SECTION 14 - TRANSPORT INFORMATION

USDOT

- No information available

Canadian TDG

- No information available.

SECTION 15 - REGULATORY INFORMATION**US FEDERAL**

TSCA

- CAS# 7732-18-5 is listed on the TSCA inventory.
- CAS# 7773-01-5 is listed on the TSCA inventory.

Health & Safety Reporting List

- None of the chemicals are on the Health & Safety Reporting List.

Chemical Test Rules

- None of the chemicals in this product are under a Chemical Test Rule.

Section 12b

- None of the chemicals are listed under TSCA Section 12b.

TSCA Significant New Use Rule

- None of the chemicals in this material have a SNUR under TSCA.

SARA

CERCLA Hazardous Substances and corresponding RQs

- None of the chemicals in this material have an RQ.

SARA Section 302 Extremely Hazardous Substances

- None of the chemicals in this product have a TPO.

Section 313

- This material contains Manganese chloride (listed as ** undefined **), 12.5%, (CAS# 7773-01-5) which is subject to the reporting requirements of Section 313 of SARA Title III and 40 CFR Part 372.

Clean Air Act:

- CAS# 7773-01-5 listed as ** no name ** is listed as a hazardous air pollutant (HAP).
- This material does not contain any Class 1 Ozone depleters.
- This material does not contain any Class 2 Ozone depleters.

Clean Water Act:

- None of the chemicals in this product are listed as Hazardous Substances under the CWA.
- None of the chemicals in this product are listed as Priority Pollutants under the CWA.
- None of the chemicals in this product are listed as Toxic Pollutants under the CWA.

OSHA:

- None of the chemicals in this product are considered highly hazardous by OSHA.

STATE

Water is not present on state lists from CA, PA, MN, MA, FL, or NJ.

Manganese chloride can be found on the following state right to know lists: California, (listed as ** no name **), Pennsylvania, (listed as ** no name **), Minnesota, (listed as ** no name **).

California No Significant Risk Level: None of the chemicals in this product are listed.

European/International Regulations

European Labeling in Accordance with EC Directives

- Hazard Symbols: XN
- Risk Phrases:
 - R 22 Harmful if swallowed.

- Safety Phrases:

WGK (Water Danger/Protection)

- CAS# 7732-18-5: No information available.
- CAS# 7773-01-5: 1

United Kingdom Occupational Exposure Limits

- CAS# 7773-01-5: OES-United Kingdom, TWA (listed as ** undefined **): 5 mg/m³ TWA (as Mn)

United Kingdom Maximum Exposure Limits

- CAS# 7773-01-5: MEL-United Kingdom, TWA (listed as ** undefined **): 0.5 mg/m³ TWA

Canada

- CAS# 7732-18-5 is listed on Canada's DSL List.
- CAS# 7773-01-5 is listed on Canada's DSL List
- This product has a WHMIS classification of D2B.
- CAS# 7732-18-5 is not listed on Canada's Ingredient Disclosure List.
- CAS# 7773-01-5 is listed on Canada's Ingredient Disclosure List.

Exposure Limits

- CAS# 7773-01-5: OEL-AUSTRALIA:TWA 5 mg(Mn)/m³ JANUARY 1993
- OEL-BELGIUM:TWA 5 mg(Mn)/m³ JANUARY 1993
- OEL-CZECHOSLOVAKIA:TWA 2 mg(Mn)/m³;STEL 6 mg(Mn)/m³ JANUARY 1993
- OEL-DENMARK:TWA 2.5 mg(Mn)/m³ JANUARY 1993
- OEL-FINLAND:TWA 2.5 mg(Mn)/m³ JANUARY 1993
- OEL-HUNGARY:TWA 0.3 mg(Mn)/m³;STEL 0.6 mg(Mn)/m³ JANUARY 1993
- OEL-JAPAN:TWA 0.3 mg(Mn)/m³ JANUARY 1993
- OEL-THE NETHERLANDS:TWA 1 mg(Mn)/m³ JANUARY 1993
- OEL-POLAND:TWA 0.3 mg(Mn)/m³ JANUARY 1993
- OEL-SWEDEN:TWA 1 mg(Mn)/m³;STEL 2.5 mg(Mn)/m³ (resp. dust)
- OEL-SWEDEN:TWA 2.5 mg(Mn)/m³;STEL 5 mg(Mn)/m³ (total dust)
- OEL-UNITED KINGDOM:TWA 5 mg(Mn)/m³ JANUARY 1993 OEL IN BULGARIA, COLOMBIA, JORDAN, KOREA check ACGIH TLV OEL IN NEW ZEALAND, SINGAPORE, VIETNAM check ACGI TLV

SECTION 16 - ADDITIONAL INFORMATION

MSDS Creation Date: 6/30/1999, Revision #4 Date: 9/02/2004

The information above is believed to be accurate and represents the best information currently available to us. However, we make no warranty of merchantability or any other warranty, express or implied, with respect to such information, and we assume no liability resulting from its use. Users should make their own investigations to determine the suitability of the information for their particular purposes. In no way shall the company be liable for any claims, losses, or damages of any third party or for lost profits or any special, indirect, incidental, consequential or exemplary damages, howsoever arising, even if the company has been advised of the possibility of such damages.

Material Safety Data Sheet

Section 1. Product and Company Identification

Product Name	Sodium Chloride	Product Code	XX0778
Manufacturer	EMD Chemicals Inc. P.O. Box 70 480 Democrat Road Gibbstown, NJ 08027 Prior to January 1, 2003 EMD Chemicals Inc. was EM Industries, Inc. or EM Science, Division of EM Industries, Inc.	Effective Date	3/3/2003
For More Information Call	856-423-6300 Technical Service Monday-Friday: 8:00 AM - 5:00 PM	In Case of Emergency Call	800-424-9300 CHEMTREC (USA) 613-996-6666 CANUTEC (Canada) 24 Hours/Day: 7 Days/Week
Synonym	SALT, HALITE		
Material Uses	Analytical reagent.		
Chemical Family	Inorganic salt.		

Section 2. Composition and Information on Ingredients

Component	CAS #	% by Weight
SODIUM CHLORIDE	7647-14-5	100

Section 3. Hazards Identification

Physical State and Appearance Solid. (Granular solid. Crystals solid.)

Emergency Overview CAUTION!
MAY CAUSE EYE AND SKIN IRRITATION.
MAY CAUSE DAMAGE TO THE FOLLOWING ORGANS: SKIN, EYES, STOMACH.

Routes of Entry Inhalation. Ingestion.

Potential Acute Health Effects
Eyes May be hazardous in case of eye contact (irritant).
Skin May be hazardous in case of skin contact (irritant). Skin inflammation is characterized by itching, scaling, reddening, or, occasionally, blistering.
Inhalation Non-hazardous in case of inhalation.
Ingestion Non-hazardous in the case of ingestion.

Potential Chronic Health Effects
Carcinogenic Effects This material is not known to cause cancer in animals or humans.

Medical Conditions Aggravated by Overexposure: Additional information See Toxicological Information (section 11)
Repeated or prolonged exposure is not known to aggravate medical condition.

Section 4. First Aid Measures

Eye Contact Check for and remove any contact lenses. In case of contact, immediately flush eyes with plenty of water for at least 15 minutes. Cold water may be used. Get medical attention.

Skin Contact In case of contact, immediately flush skin with plenty of water. Cover the irritated skin with an emollient. Remove contaminated clothing and shoes. Cold water may be used. Wash clothing before reuse. Thoroughly clean shoes before reuse. Get medical attention.

Inhalation If inhaled, remove to fresh air. If not breathing, give artificial respiration. If breathing is difficult, give oxygen. Get medical attention.

Ingestion Do NOT induce vomiting unless directed to do so by medical personnel. Never give anything by mouth to an unconscious person. If large quantities of this material are swallowed, call a physician immediately. Loosen tight clothing such as a collar, tie, belt or waistband.

Section 5. Fire Fighting Measures

Flammability of the Product Non-flammable.

Auto-ignition Temperature Not applicable.

Flash Points Not applicable.

Flammable Limits Not applicable.

Products of Combustion Not available.

Fire Hazards in Presence of Various Substances Not applicable.

Explosion Hazards in Presence of Various Substances Risks of explosion of the product in presence of static discharge: No.

Fire Fighting Media and Instructions Risks of explosion of the product in presence of mechanical impact: No.
Not applicable.

Protective Clothing (Fire Hazards) Not applicable.

Special Remarks on Fire Hazards Not available.

Special Remarks on Explosion Hazards Not available.

Section 6. Accidental Release Measures

Small Spill and Leak	Use appropriate tools to put the spilled solid in a convenient waste disposal container.
Large Spill and Leak	Use a shovel to put the material into a convenient waste disposal container. Finish cleaning by spreading water on the contaminated surface and allow to evacuate through the sanitary system.
Spill Kit Information	No specific spill kit required for this product.

Section 7. Handling and Storage

Handling	Avoid contact with eyes, skin and clothing. Do not ingest.
Storage	Keep container tightly closed. Keep container in a cool, well-ventilated area.

Section 8. Exposure Controls/Personal Protection

Engineering Controls Use process enclosures, local exhaust ventilation, or other engineering controls to keep airborne levels below recommended exposure limits. If user operations generate dust, fume or mist, use ventilation to keep exposure to airborne contaminants below the exposure limit.

Personal Protection	Eyes Splash goggles. Body Lab coat. Respiratory Dust respirator. Be sure to use an approved/certified respirator or equivalent. Hands Gloves. Feet Not applicable.
Protective Clothing (Pictograms)	
Personal Protection in Case of a Large Spill	Splash goggles. Full suit. Dust respirator. Boots. Gloves. A self-contained breathing apparatus should be used to avoid inhalation of the product. Suggested protective clothing might not be sufficient; consult a specialist BEFORE handling this product.
Product Name	Exposure Limits
SODIUM CHLORIDE	Not available.

Section 9. Physical and Chemical Properties

Odor	Odorless.
Color	White.
Physical State and Appearance	Solid. (Granular solid. Crystals solid.)
Molecular Weight	58.44 g/mole
Molecular Formula	ClNa
pH	Not available.
Boiling/Condensation Point	1430.9°C (2607.6°F)
Melting/Freezing Point	800.9°C (1473.6°F)
Specific Gravity	Not available.
Vapor Pressure	Not available.
Vapor Density	Not available.
Odor Threshold	Not available.
Evaporation Rate	Not available.
LogKow	Not available.
Solubility	Soluble in water.

Section 10. Stability and Reactivity

Stability and Reactivity	The product is stable.
Conditions of Instability	Not available.
Incompatibility with Various Substances	Not available.
Reactivity with Water	Not available.
Hazardous Decomposition Products	Not available.
Hazardous Polymerization	Will not occur.

Section 11. Toxicological Information

RTECS Number:	Sodium Chloride	VZ4725000
Toxicity	Acute oral toxicity (LD50): 3000 mg/kg [Rat]. Acute toxicity of the vapor (LC50): >42000 mg/m ³ 1 hour(s) [Rat].	
Chronic Effects on Humans	Contains material which may cause damage to the following organs: skin, eyes, stomach.	
Acute Effects on Humans	May be hazardous in case of eye contact (irritant). May be hazardous in case of skin contact (irritant). Skin inflammation is characterized by itching, scaling, reddening, or, occasionally, blistering. Non-hazardous in case of inhalation. Non-hazardous in the case of ingestion.	
Synergistic Products (Toxicologically)	Not available.	
Irritancy	Draize Test (Rabbit): Eye: 100mg/24h. Reaction: Moderate. Skin: 500 mg/24h mild	
Sensitization	Not available.	
Carcinogenic Effects	This material is not known to cause cancer in animals or humans.	
Toxicity to Reproductive System	Tests on laboratory animals for reproductive effects are cited in Registry of Toxic Effects on Chemical Substances (RTECS).	
Teratogenic Effects	Not available.	
Mutagenic Effects	Tests on laboratory animals for mutagenic effects are cited in Registry of Toxic Effects of Chemical Substances	

(RTECS).

Section 12. Ecological Information

Ecotoxicity Not available.
 BOD5 and COD Not available.
 Toxicity of the Products of Biodegradation The product itself and its products of degradation are not toxic.

Section 13. Disposal Considerations

EPA Waste Number Not available.
 Treatment Material does not have an EPA Waste Number and is not a listed waste, however consultation with a permitted waste disposal site (TSD) should be accomplished. Always contact a permitted waste disposal (TSD) to assure compliance with all current local, state, and Federal Regulations.

Section 14. Transport Information

DOT Classification Not available.
 TDG Classification Not available.
 IMO/MDG Classification Not available.
 ICAO/IATA Classification Not available.

Section 15. Regulatory Information

U.S. Federal Regulations TSCA 8(b) inventory: SODIUM CHLORIDE
 SARA 302/304/311/312 extremely hazardous substances: No products were found.
 SARA 302/304 emergency planning and notification: No products were found.
 SARA 302/304/311/312 hazardous chemicals: SODIUM CHLORIDE
 SARA 311/312 MSDS distribution - chemical inventory - hazard identification: SODIUM CHLORIDE:
 Immediate (Acute) Health Hazard, Delayed (Chronic) Health Hazard
 SARA 313 toxic chemical notification and release reporting: No products were found.
 Clean Water Act (CWA) 307: No products were found.
 Clean Water Act (CWA) 311: No products were found.
 Clean air act (CAA) 112 accidental release prevention: No products were found.
 Clean air act (CAA) 112 regulated flammable substances: No products were found.
 Clean air act (CAA) 112 regulated toxic substances: No products were found.
 WHMIS (Canada) Not controlled under WHMIS (Canada).
 CEPA DSL: SODIUM CHLORIDE
 This product has been classified in accordance with the hazard criteria of the Controlled Product Regulations and the MSDS contains all required information.
 International Regulations
 EINECS SODIUM CHLORIDE 231-598-3
 DSCL (EEC) R36/38- Irritating to eyes and skin.
 International Lists Australia (NICNAS): SODIUM CHLORIDE
 Japan (MITI): SODIUM CHLORIDE
 Korea (TCCL): SODIUM CHLORIDE
 Philippines (RA6969): SODIUM CHLORIDE
 China: No products were found.
 State Regulations No products were found.
 California prop. 65: No products were found.

Section 16. Other Information



National Fire Protection Association (U.S.A.)
 Fire Hazard 0
 Health 0 0
 Reactivity
 Specific Hazard


Changed Since Last Revision +


Notice to Reader

The statements contained herein are based upon technical data that EMD Chemicals Inc. believes to be reliable, are offered for information purposes only and as a guide to the appropriate precautionary and emergency handling of the material by a properly trained person having the necessary technical skills. Users should consider these data only as a supplement to other information gathered by them and must make independent determinations of suitability and completeness of information from all sources to assure proper use, storage and disposal of these materials and the safety and health of employees and customers and the protection of the environment. EMD CHEMICALS INC. MAKES NO REPRESENTATION OR WARRANTY OF ANY KIND, EXPRESS OR IMPLIED, INCLUDING MERCHANTABILITY OR FITNESS FOR A PARTICULAR USE, WITH RESPECT TO THE INFORMATION HEREIN OR THE PRODUCT TO WHICH THE INFORMATION REFERS.

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
Artéfact: ringing artifact**Severe ringing**

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Artéfact: sagittal ghosts**Sagittal ghosts**

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